

Effect of Olive Oil Adulteration on Peroxide Value, Delta-K and on the Acidity Nabali-Baladi Olive Oil Quality

Orwa Jaber Houshia^{1,*}, Oday Zaid², Hazem Shqair², Motasem Zaid², Nashida Fashafsheh¹,
Reeman Bzoor¹

¹Arab American University, P.O. Box 240, Jenin, Palestine
²National Agricultural Research Center (NARC), Jenin, Palestine

Abstract The purpose of the study was to evaluate the effects of olive oil adulteration on the peroxide value and oil acidity, and also to evaluate the effects of the olive oil adulteration on the delta K value, which classifies the purity (grade) of olive oil. The study was begun by collecting several samples of olive oil from northern governorates of Palestine. Also, sunflower, corn, and soybean oil samples were purchased at the local market. Olive oil samples were analyzed for delta K value, peroxide value and acidity. Olive oil samples were analyzed at 232, 268, 270 and 274 nm using UV-VIS instrument. Then, the samples were mixed with other oils at certain ratios and evaluated at the same wavelength. All samples were analyzed using quartz cuvettes. Cyclohexane was used as solvent for the oil. The same samples were evaluated in the same manner by titrating the samples with thiosulfate to get the peroxide value, and also with potassium hydroxide to get the acidity of oils. The IOC standard is <20 mEq O₂ / kg oil. The acidity test (the IOC standard for free fatty acids in extra virgin olive oil is a maximum of 0.8%.) can differentiate between authentic versus adulterated oil. Reliable technique was used for measuring the authenticity of olive oil.

Keywords Nabali Baladi, Olive Oil Adulteration, Peroxide Value, Acidity, Delta K

1. Introduction

Olive oil is one of the oldest known vegetable oils mainly produced in the countries surrounding the Mediterranean Sea. It is a natural fruit juice, obtained from the fruit of the Olive tree, with a unique composition and quality. Olive oil is an integral part of the Palestinian economy. About 95% of the olives picked in the Palestinian territories is converted to olive oil, most of which is exported internationally. Olive oil plays an essential part of the Palestinian and the Mediterranean diet. This oil has several complex chemical components. The composition of olive oil can give valuable information in understanding their functional, quality and nutritional properties. The main type of fat found in all kinds of olive oil is monounsaturated fatty acids (MUFAs). MUFAs are actually considered a healthy dietary fat. If your diet emphasizes unsaturated fats, such as MUFAs and polyunsaturated fats (PUFAs), instead of saturated fats and Trans fats, you may gain certain health benefits. Olive oil is one of the very few oils that can be consumed in its natural

form, thus preserving all its natural constituents [1, 2].

Olive oil is well established in the food industry. Demand continues to grow not only because of its distinct flavor, but also because of an increased awareness of its health benefits. The increase in the world price of olive oil compared to other vegetable oils make it Susceptible to fraud commonly adulteration with cheap vegetable oils [3]. For this, the International Olive oil Council (IOC) placed strict criteria for olive oil classifications by chemical, sensory and organoleptic properties [4-6]. Thus, Olive oil types [7] can be also classified as:

- Extra virgin olive oil: virgin olive oil with a free acidity, expressed as oleic acid, of not more than 0.8 grams per 100 grams.
- Virgin olive oil: virgin olive oil with a free acidity, expressed as oleic acid, of not more than 2.0 grams per 100 grams.
- Ordinary virgin olive oil: virgin olive oil with a free acidity, expressed as oleic acid, of not more than 3.3 grams per 100 grams.
- Refined olive oil: olive oil obtained from virgin olive oils by refining methods which do not lead to alterations in the initial glyceridic structure. It has a free acidity, expressed as oleic acid, of not more than 0.3 grams per 100 grams.

* Corresponding author:

orwa.houshia@aauj.edu (Orwa Jaber Houshia)

Published online at <http://journal.sapub.org/als>

Copyright © 2014 Scientific & Academic Publishing. All Rights Reserved

There are several methods of analysis that have been proposed and are used to evaluate the unique olive oil composition and quality. The detection of olive oil adulteration is often complicated with no single test that can accomplish the task. A series of tests are employed to determine Olive oil authenticity, quality and identity of the adulterant. Included in this testing method is the determination of free acidity, peroxide value, UV extinction, fatty acid composition, and the Bellier test. Test results are measured against the International Olive Council trade standard to identify abnormalities. Each test provides key information which allows a decision to be made with respect to the grade of Olive oil and the identity of any adulterants [8].

Previous researchers were able distinguish adulterated olive oil from extra virgin olive by using ^{31}P NMR [9], however Spectrophotometric examination in the ultraviolet can provide information on the quality of a fat. The absorption at the wavelengths specified in the method is due to the presence of conjugated diene and triene systems. These absorptions are expressed as specific extinctions $E_{1\%}$ (the extinction of 1% solution of the fat in the specified solvent, in a thickness of 1 cm) conventionally indicated by K (also referred to as "extinction coefficient". A low absorption in this region is indicative of a high quality extra virgin olive oil, whereas adulterated/refined oils show a greater level of absorption in this region. Hence arose the idea of our research based on the study of changes in the values of absorption for pure virgin olive oil with a certain mixed percentage of vegetable oils at 232 nm spectrophotometrically. And in the ultraviolet domain nm, 270 nm (270nm rate of change in the values of absorption a) and to determine ΔK value by measuring the absorbance at 266nm, and the 274 nm in order to determine the lowest rates of cheating vegetable oil until it is exposed.

2. Methodology

Three samples were purchased for each of the extra virgin olive oil, sunflower oil, corn oil, Soybean oil and cottonseed oil randomly from the local markets and different sources, taking into account the availability of brands and the shelf life on the package. Pure olive oil were from Jalboon and Kefret areas from the westbank.

Adulterated samples of pure olive oil mixed with corn, soy, sunflower oil, were made at these percentages: 5%, 10%, 20%, 30%, and 40%. And a sample of pure olive oil 0% adulteration to serve as a reference. Samples well-kept in a dark place and in airtight containers until the end of the analysis.

Delta (Δ) k is evaluated using the following formula:

$$\Delta k = \frac{k_{266} + k_{274}}{2} - k_{270}$$

Where:

E : absorbance
 $k = \frac{E}{c.s.}$: concentration
 s : cuvette path length
 $R\text{-value} = \frac{k_{232\text{nm}}}{k_{270\text{nm}}}$

2.1. Spectrophotometric Measurements

Perkin-Elmer spectrophotometer with ability measuring extinction in the ultraviolet between 220 and 360nm. Rectangular quartz cuvettes, with covers, having an optical length of 1 cm were used. The method is most often used in a quantitative way to determine concentrations of an absorbing species in solution, using the Beer-Lambert law

Olive oil samples (50 in total) were diluted in cyclohexane. All samples were measured in matched, synthetic fused silica cuvettes (10 mm is the recommended path length). Weigh accurately 0.25 g of the sample so prepared into a 25 ml graduated flask, make up to the mark with the solvent specified and homogenize. The resulting solution must be perfectly clear. If opalescence or turbidity is present, filter quickly through paper. Fill a cuvette with the solution obtained and measure the extinctions at an appropriate wavelength a 232, 266, 270 and 274 nm, using the solvent used as a reference.

2.2. Peroxide Value

Quantity of those substances in the sample, expressed in terms of active oxygen, which oxidize potassium iodide under the conditions specified in this International Standard, divided by the mass of the test portion. Add 50 ml of acetic acid/isooctane solution to 10g of oil in conical flask and replace the stopper. Swirl the flask until the sample has dissolved. Using a suitable volumetric pipette, add 0.5 ml of saturated potassium iodide solution and replace the stopper. Allow the solution to react for 1 min, thoroughly shaking it at least three times during this period, then immediately add 30 ml of distilled water.

Titrate the solution with sodium thiosulfate solution, adding it gradually and with constant, vigorous agitation, until the yellow iodine color has almost disappeared. Add about 0.5 ml of starch solution and continue the titration with constant agitation, especially near the endpoint, to liberate all of the iodine from the solvent layer, adding the sodium thiosulfate solution drop by drop until the blue color just disappears.

2.3. Acidity Value

Weigh into a flask a sufficient mass of the test sample as shown in **Table 1**, according to the color and expected acid value. NOTE: The mass of the test portion and the concentration of the titrate should be such that the titrate does not exceed 10 ml.

Heat to boiling 50 ml of the ethanol containing 0.5 ml of the phenolphthalein indicator in a second flask. Whilst the temperature of the ethanol is still over 70°C, neutralize it carefully with a solution of 0.1 mol/l sodium or potassium hydroxide. The endpoint of the titration is reached when the


addition of a single drop of alkali produces a slight but definite color change persisting for at least 15 s. Add the neutralized ethanol to the test portion in the first flask and mix thoroughly. Bring the contents to the boil and titrate with the sodium or potassium hydroxide solution, (depending on the expected acidity of the sample), agitating the flask contents vigorously during the titration.

Table 1. Mass of test portion and accuracy of weighing

Expected acid value	Mass in (g)	Weighing Accuracy (g)
< 1	20	0.05
1-4	10	0.02
4-15	2.5	0.01

3. Results and Discussion

The quality of the olive oil is studied by measuring the characteristics of the absorption bands between 200 and 300 nm. These are frequencies related to conjugated diene and triene systems. A low absorption in this region is indicative of a high quality extra virgin olive oil, whereas adulterated/refined oils show a greater level of absorption in this region. The absorbency at 232 nm is caused by hydroperoxides (primary stage of oxidation) and conjugated dienes (intermediate stage of oxidation). The absorbency at 270 nm is caused by carbonylic compounds (secondary stage of oxidation) and conjugated trienes (technological treatments). The index “(Δ) k” is a criterion of discrimination between a bad quality virgin olive oil and a virgin olive oil adulterated with refined olive oil. The degree of oxidation of olive oil is reflected by its specific extinction at 232 nm and 270 nm. In fact, virgin olive oils, of good quality and correctly stored, contain very few products of oxidation; these mainly peroxidic in nature, have a maximum absorption at approximately 232 nm. The values of “(Δ) k”, at 232 and 270 nm in such olive oils are below the maximum provided for in the standard. On the other hand, when the oil is treated with a decolorizing agent (i.e. an absorbent earth) during the refining process, conjugated trienoic compounds are formed. These compounds have a maximum absorption situated at approximately 270 nm; this means that refined oils have higher values of “k” at 270 nm. The absorption in the U.V. is related to the presence of conjugate double bonds. In the oils, due to oxygen fixation in linolenic and linoleic acids’ double bond position, hydroperoxides arise. The double bond nature that we can obtain provoke the formation of conjugate diene systems between carbon atoms.



This kind of conjugate systems presents a maximum absorption at 232 nm wavelength. Instead, during more advanced oxidation states, products are generated with conjugate diene systems carbon-oxygen like:



The maximum absorption in this case ranges between wavelength of 260/280 nm. A conjugate triene system presents a triple band with a maximum absorption at 268 nm and another at 232 nm. The absorption curve in the U.V. of oil is very much influenced by the "oxidation" products, some of which provoke an increase of the absorption at 232 nm and others at 270 nm. This curve is called Delta (Δ) k and is part of the fundamental analysis that we conducted to establish the purity and degree of degradation of an oil Delta (Δ) k correlates with the state of oxidation by detecting specific oxidized compounds, some generated from secondary oxidation, and also detect possible adulteration with refined oils. Oxidation produces peroxides from some fatty acids. Autoxidation (when the oil gets in contact with oxygen in the air, first during the milling process, and later on during storage in tanks and bottles) is the most common process leading to oxidative deterioration and is defined as the spontaneous reaction of atmospheric oxygen with lipids. [10]

Peroxides are unstable compounds which are further oxidized to yield volatile and non-volatile components that give rise to off-flavors and undesirable aromas in the oil. This is referred to as secondary stage oxidation. Oxidation is accelerated at higher temperatures, which is called thermal oxidation, with increases in free fatty acid and polar matter contents, foaming, color, and viscosity [11]. The autoxidation and thermal oxidation of unsaturated fatty acids occurs via a free radical chain reaction that proceeds through three steps of initiation, propagation, and termination. A simplified scheme explaining the mechanism of autoxidation is shown in **Figure 1**.

As oxidation normally proceeds very slowly at the initial stage, the time to reach a sudden increase in oxidation rate is referred to as the induction period [12]. Lipid hydroperoxides have been identified as primary products of autoxidation; decomposition of hydroperoxides yields aldehydes, ketones, alcohols, hydrocarbons, volatile organic acids, and epoxy compounds, known as secondary oxidation products. These compounds, together with free radicals, constitute the bases for measurement of oxidative deterioration of food lipids.

The peroxide value (PV) is an indicator of the initial stages of oxidative change [13]. The PV represents the total hydroperoxide content and is one of the most common quality indicators of fats and oils during production and storage [14, 15]. A number of methods have been developed for determination of PV, among which the iodometric titration, ferric ion complex measurement spectrophotometry, and infrared spectroscopy are most frequently used.

Iodometric titration assay, which is based on the oxidation of the iodide ion (I-) by hydroperoxides (ROOH), is the basis of current standard methods for determination of PV. In this method, a saturated solution of potassium iodide is added to

oil samples to react with hydroperoxides. The liberated iodine (I_2) is then titrated with a standardized solution of sodium thiosulfate and starch as an endpoint indicator. The PV is obtained by calculation and reported as milliequivalents of oxygen per kilogram of sample (meq/kg). The official determination is described by IUPAC [16].

Chemical reactions involved are shown in **Figure 2**. The assay includes six steps: accurate weighing of the sample, dissolution of lipids in chloroform, acidification with acetic acid, addition of potassium iodide, incubation for exactly 5 minutes, and titration with sodium thiosulfate.




Figure 1. Proposed mechanism of autoxidation




Figure 2. Chemical reactions involved in iodometric titration of olive oil




Figure 3. Oil hydrolysis with KOH

As mentioned earlier, free fatty acid formation precedes most, if not all oil deterioration. Free acidity levels increase by hydrolysis of the major oil molecules in the fruit (triglycerides) early in the production process, from harvest through milling, while water and plant enzymes are still in contact with the oil. For this reason, milling the olives soon after harvest, and promptly separating oil from vegetation water are critical to maintaining a low initial free acidity level, essential for higher-quality oil. Free acidity values provide an indication of how the fruit was handled prior to processing and the length of time from harvest to milling. Free acidity is also an early indicator of the potential longevity of the oil. Higher quality oils recently produced, will exhibit very low acidity—it may be no higher than 0.35% in the best extra virgin olive oils (EVOO).

The acid value is the number of milligrams of potassium hydroxide required to neutralize the free fatty acids present in 1 g of fat, when determined in accordance with the procedure specified in this International Standard. Acid value is expressed in milligrams per gram. However, acidity is defined as the content of free fatty acids determined according to the procedure specified in this International

Standard. Acidity is expressed as a percentage by mass. The hydrolysis of the major oil molecules as it reacts with potassium hydroxide is shown in **Figure 3**, it gives three molecules of free fatty acids (FFA).

Acidity level is a measure of the free fatty acids in the oil, which essentially indicates how ripe, or more accurately, how decomposed the oil is when it went into the bottle. The IOC standard for free fatty acids in extra virgin olive oil is a maximum of 0.8%, but fresh, high-quality EVOO can certainly achieve FFA values of 0.35% or lower.

The Palestinian National Authority recently set the following classification for olive quality:

1. Premium Extra Virgin: acidity < 1%, and PV < 12
2. Premium Virgin: acidity < 1.5% and PV < 20
3. Premium Pseudo-Virgin: acidity < 3.3% and PV < 20
4. Lampant Virgin: acidity > 3.3% and PV > 20

The effects of percent mixing olive oil with other oil on the peroxide value and acidity value was studied, and the selected oils was corn, sun flower and soy. The following graph shows the percent mixing on the y-axis and peroxide value on the x-axis.




Figure 4. Peroxide value of % mixing of olive oil with three type vegetable oil: PON = soy oil, POC = corn oil, POS = sunflower oil. Kferet samples

Table 2. Peroxide value for mixing olive oil from different places

	% mixing	0	5	10	20	40	100	time
PV of kferet	P.O.S	27.68	21.97	22.86	20.26	18.49	8.08	1 st week
	P.O.C	27.68	20.27	20.7	22.7	17.25	4.89	2 nd week
	P.O.N	22.78	21.88	20.4	19.54	16	3.8	2 nd week
PV of jalbone	P.O.S	22.45	22.07	21.19	20.06	15.89	5.7	1 st week
	P.O.C	27.28	22.08	24.27	46.55	12.99	7.29	1 st week
	P.O.N	27.28	23.27	25.17	22.98	21.88	3.2	2 nd week

In general, all trends of peroxide value for three % mixing oil down slop. Its observed that when the % mixing is increase the peroxide value of Kferet olive oil is decrease. from the graph we can see that the peroxide value of pure olive oil is very high above IOC standard that because it is old oil from two years but when we mix it with other oil (adulterate) the peroxide value was lowered and its seem clearly on the higher mixing% as 20, 40, and sometime 10%, but the adulteration with 5% is not different too much from 0% (pure olive oil) in most sample that we tests. The same thing for jalbon olive oil as it appear in figure 5.

The change in peroxide value with deferent %mixing is observed in next table. The table show the %mixing and peroxide value for three mixing oil sun flower(S), corn(C), and soya oil(N) for kferet and Jalbone village.

The effects of % mixing olive oil with other oil on acidity

value was studied, and the selected oils was corn, sun flower and soy. The following graph shows the % mixing on the y-axis and acidity value on the x-axis for the kferet samples.

The same thing, all trends of acidity value for three % mixing oil down slop. Its observed that when the % mixing is increase the acidity value of Kferet olive oil is decrease. from the graph in figure 6, we can see that the acidity value of pure olive oil is very high above IOC standard that because it is old oil from two years so there is a large amount of free fatty acid but when we mix it with other oil (adulterate) the acidity value was lowered and its seem clearly on the higher mixing % as 20, 40, and sometime 10%, but the adulteration with 5% is not different too much from 0% (pure olive oil) in most sample that we tests. The same thing for jalbon olive oil as it appear in figure 7.




Figure 5. Peroxide value of % mixing Jalboon olive oil with three type mixing oil




Figure 6. Acidity value of % mixing olive oil from Kferet village for three mixing oil




Figure 7. Acidity value of % mixing Jalbone olive oil with three type mixing oil

Table 3. Acidity value for mixing olive oil from different places

		% mixing	0	5	10	20	40	100	time
AV of kferet	A.O.S		3.66	3.55	3.41	2.97	2.27	0.14	1 st week
	A.O.C		3.66	5.09	3.41	2.99	2.35	0.17	2 nd week
	A.O.N		3.92	3.8	3.41	3.05	2.24	0.17	2 nd week
AV of jalbone	A.O.S		1.37	1.42	1.23	1.06	0.9	0.11	1 st week
	A.O.C		1.37	1.34	1.23	1.12	0.89	0.17	1 st week
	A.O.N		1.37	1.34	1.37	1.12	0.9	0.11	2 nd week




Figure 8. Percent adulteration verses Delta-K

The change in acidity value with deferent %mixing is observed in next table. The table show the %mixing and acidity value for three mixing oil sun flower(S), corn(C), and soya oil(N) for kferet and Jalbone village with the time of the tests work.

Figure 8 show the relationship between % of mixing Jalboon olive oil with soya oil and delta k. Its appear a positive and linear relationship between % of mixing and delta k with no deviation, when % of mixing increase the

delta k value increase.

The absorbance at the four λ s of all samples is presented at y axis and percent of adulteration is presented at x axis, so three graphs are obtained as shown in figures 9, 10 and 11

It is appear form these graphs in figures 9, 10, and 11 that the absorbance at λ 232nm at all percent of mixing is nearly the same, but at 266 nm, 270 nm and 274 nm the absorbance is increase as the percent of mixing increase except at 5% is decrease and it the lowest point.




Figure 9. Percent olive oil mixing with sun-flower oil verses change in absorbance




Figure 10. Percent olive oil mixing with soy-bean oil verses change in absorbance




Figure 11. Percent olive oil mixing with corn oil verses change in absorbance

From these results at λ 232 can't be used to detect the oil adulteration because the absorbance at all percent of mixing is nearly the same, but λ 266, λ 270 and λ 274 can be used to detect oil adulteration because the absorbance increase when the percent of adult increase but above 5%.

So this method is not effective at 5% so we need another way to detect adulteration at 5% and less.

4. Conclusions

The index Δk is very important to distinguish between bad quality origin olive oil and origin olive oil adulterated with refined olive oil.

From the results, when the percent of mixing olive oil with other oils increase, delta k increase, and the absorbance at 232 nm at all percent of mixing is nearly the same with slightly difference so it can't be used to know if the oil is adult or not, but at the other percents the absorbance increase as the percent increase except at 5% it decrease so these wavelength can be used to know if the oil is adult or not but above 5%.

It appears that mixing olive oil with other oil had no significant effect to improve these tests to detect the adulteration, because at small % mixing the value of peroxide and acidity value did not give a clear result about the adulteration or quality of olive oil because it is near from the 0% result, and at higher % mixing it didn't give the real quality of oil because as mixing % increase peroxide and acidity value decrease, and that what is told by curves, which show that mixing olive oils actually improves peroxide and acidity, so we can't depend on these tests by 100% to judge the quality of olive oil.

ACKNOWLEDGEMENTS

The authors would like to thank the AAUJ chemistry undergraduate students and NARC scientists for their collaboration and team work.

REFERENCES

- [1] Zaid, O. Houshia, O.J., AbuEid, M., Zaid, M., 2013. Palestinian Nabali-Baladi Olive Oil Quality: Premium Ultra Fine Extra Virgin Olive Oil Classification. USARJ, Vol. 1, No. 2, pp. 29-34,
- [2] Aparicio, R. and Aparicio-ruiz, R. (2000). Authentication of Vegetable Oils by Chromatographic Techniques. *J. Chromatography A*, 881(0021-9673), 93-104.
- [3] Hidalgo, C. F., Navas, M. A. and Garin, A. (1993). Virgin Olive Oil Quality, Some New Criteria for its Evaluation. *Grasas-y-Aceites*. (Espana). 44(1), 10-17.
- [4] Gonzalez, D. R. and Aparicio, R. (2002). Detection of Defective Virgin Oils by Metal Oxide Sensors. *European Food Research and Technology*. Instituto de la Grasa, Spain.
- [5] Oliveros, C. M., Pavon, J. L., Pinto, C. G., Laespada, F. M. E., Cordero, B. M. and Forina, M. (2002). Electronic Nose Based on Metal Oxide Semiconductor Sensors as a Fast Alternative for the Detection of Adulteration of Virgin Olive Oils. *Universidad de Salamanca, Spain, Università degli Studi di Genova, Italy Analytica Chimica Acta*, 459 (2,22), 219-228.
- [6] INTERNATIONAL OLIVE COUNCIL, COI/T20/Doc. no. 19/Rev. 3, 2010, METHOD OF ANALYSIS SPECTROPHOTOMETRIC INVESTIGATION IN THE ULTRAVIOLET <http://www.internationaloliveoil.org/>.
- [7] Z. Piravi Vanak, Jahan B. Ghasemi, M. Ghavami and H. Ezzatpanah 2010 "Detection and Quantification of Adulteration in Olive Oils by Global Method and Extinction Coefficient", *Australian Journal of Basic and Applied Sciences*, 4(12): 6254 6259.
- [8] Determination of Peroxide Value, in IUPAC: Standard Methods for the Analysis of Oils, Fats and Derivatives, IUPAC Standard Method 2.501, International Union of Pure and Applied Chemistry, Pergamon, Oxford, England, 1992
- [9] Georgia, F., Apostolos, S., George, S., Emmanuel, S. and Photis, D. (2005). Detection of Extra Virgin Olive Oil Adulteration with Lampante Olive Oil and Refined Olive Oil Using Nuclear Magnetic Resonance Spectroscopy and Multivariate Statistical Analysis. *J. Agric. Food Chem.* 53, 2810-2816.
- [10] M. Gordon, in J. Pokorny, N. Yanishlieva, and M. Gordon, eds., *Antioxidants in Food: Practical Applications*, Woodhead Publishing, Ltd., Cambridge, England, 2001, pp. 7-21.
- [11] E. G. Perkins, in A. J. St. Angelo, ed., *Lipid Oxidation in Food*, American Chemical Society, Washington, D.C., 1992, pp. 310-321.
- [12] J. Velasco, M. L. Anderson, and L. H. Skibsted, 2004. *Food Chem.*, 85, 623-632.
- [13] J. Velasco and C. Dobarganes, *Eur. J. Lipid Sci. Technol.*, 104, 661-676 (2002).
- [14] M. Antolovich, P. D. Prenzler, E. Patsalides, S. McDonald, and K. Robards, *Analyst*, 127, 183-198 (2002).
- [15] A. Riuz, M. J. Ayora-Canada, and B. Lendl, *Analyst*, 126, 242-246 (2001).
- [16] Determination of Peroxide Value, in IUPAC: Standard Methods for the Analysis of Oils, Fats and Derivatives, IUPAC Standard Method 2.501, International Union of Pure and Applied Chemistry, Pergamon, Oxford, England, 1992.
- [17] Goodacre, R., D.B Kell, and G. Bianchi, Rapid Assessment of the Adulteration of Virgin Olive Oils by Other Seed Oils Using Pyrolysis Mass Spectrometry and Artificial Neural Networks, *J.Sci. Food Agric.* 63:297-307 (1993).
- [18] Clodoveo, M. L., Delcuratolo, D., Gomes, T., Colelli, G., Effect of different temperatures and storage atmospheres on Coratina olive oil quality. *Food Chem.* 2007, 102, 571-576.
- [19] Galtier O., N. Dupuy, Le Y. Dr eau, D. Ollivier, C. Pinatel, J.Kister and J. Artaud, 2007. Geographic origins and compositions of virgin olive oil determination by chemometric analysis of NIR spectra. *Analytica Chimica*

Acta., 595:163-144.

[20] F. R. Van de Voort, A. A. Ismail, J. Sedman, J. Dubois, T.

Nicodemo: (1994) The determination of peroxide value by
furrier transform infrared spectroscopy. J. Am. Oil Chem.
Soc. 71 921–926.